


VERSATILITY IS VITAL.

# Vital™ Deformity

Spinal Fixation System



HIGHRIDGE

# Key Instruments

## Vital Derotation

- Versatile and efficient instrumentation
- Deformity correction in all three anatomical planes




## Rod Reducers

- Multiple reducer options
- Up to 30 mm of reduction and 15 mm of medial-lateral translation

# New Offerings

## Vitalign™ S-rods

- Designed to aid sagittal alignment\*
- Save time on rod bending



Available in Ti and CoCr,  $\varnothing 5.5/6.0 \times 510$  mm  
Low, Average, and High patient PI options

- 36° of distal lumbar lordosis
- Up to 30° of proximal lumbar lordosis
- Up to 45° of thoracic kyphosis

## Rod Connectors

- Large offset to avoid existing hardware
- Added connections to interface with side hook holders



For more information, visit [highridgemedical.com](https://highridgemedical.com)



**HIGHRIDGE**

CONFIDENCE THAT MOVES

\*Data on file.

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of Highridge Medical, LLC or an affiliate; and all products are manufactured by one or more of the spinal subsidiaries of Highridge Medical, LLC and marketed and distributed by its authorized marketing partners. Please refer to the Instructions for Use, package labeling, and applicable Surgical Technique Guide. Product clearance and availability may be limited to certain countries/regions. This material does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Highridge Medical. HM0044 REV A US-EN 08/24 ©2024 Highridge Medical, LLC.